

PLEASE HANDLE WITH CARE

I am

'S wheelchair

Phone: +34



1.

Please treat me with care.



2.

Do not place anything on top of me.



3.

I do not fold and should not be dismantled.



4.

I am essential for my owner's mobility.


TECHNICAL INFORMATION



Model:



Dimensions:



Battery:



Transport certification:




HANDLING INSTRUCTIONS

- 1 The wheelchair is switched off by turning off the joystick. Do not disconnect any cables or remove the batteries.
- 2 It can be set to **MANUAL** mode or **ELECTRIC** mode (to keep it braked).
- 3 Turn the **red levers** located under the seat, one on each side.

MANUAL (DISCONNECT)


Fig. 52



CONDUCCIÓN

ELECTRIC (BRAKE)

Fig. 53



PALANCA



MANUAL mode: wheelchair can be pushed.

ELECTRIC mode: wheelchair stays braked.



BATTERY AIR TRANSPORT CERTIFICATE

CERTIFICATE OF AIR TRANSPORTATION



THIS IS TO CERTIFY THAT BATTERIES MANUFACTURED BY HAZE UNDER THE BRAND NAME **HAZE** ARE SUITABLE AND APPROVED FOR TRANSPORTATION BY ROAD, SEA AND AIR FREIGHT.

Type: AGM and GELLED ELECTROLYTE - Ranges

Voltage: 2, 4, 6 & 12 Volt

UNDER NATIONAL AND INTERNATIONAL STANDARDS SEALED MAINTENANCE FREE LEAD ACID BATTERIES ARE CLASSIFIED AS NON-DANGEROUS GOODS.

CURRENT APPLICABLE REGULATIONS

ROAD FREIGHT: US DOT 49CFR173.159(d) and ICAO/IATA special provisions.

SEA FREIGHT: IMDG, U.N.-Nr. 2800, CLASSIFICATION B

AIR FREIGHT: IATA, UN.- Nr. 2800, CLASSIFICATION 8

Authorised Signature

Thank you for taking care of me.